



US006348070B1

(12) **United States Patent**
Teissl et al.

(10) Patent No.: **US 6,348,070 B1**

(45) Date of Patent: **Feb. 19, 2002**

(54) **MAGNETIC-INTERFERENCE-FREE
SURGICAL PROSTHESES**

(75) Inventors: **Christian Teissl, Innsbruck; Erwin
Hochmair, Axams, both of (AT)**

(73) Assignee: **Med-El Elektromedizinische Gerate
Ges.m.b.H (AT)**

(*) Notice: Subject to any disclaimer, the term of this
patent is extended or adjusted under 35
U.S.C. 154(b) by 0 days.

(21) Appl. No.: **09/295,459**

(22) Filed: **Apr. 16, 1999**

Related U.S. Application Data

(60) Provisional application No. 60/082,133, filed on Apr. 17,
1998.

(51) Int. Cl.⁷ **A61F 2/02; A61F 2/18**

(52) U.S. Cl. **623/11.11; 623/10; 623/24;
600/12; 607/60**

(58) Field of Search **623/10, 24, 11.11;
600/302, 423, 424, 12; 607/60**

(56) **References Cited**

U.S. PATENT DOCUMENTS

4,038,990 A	*	8/1977	Thompson	607/9
RE32,947 E		6/1989	Dormer et al.	
4,918,745 A		4/1990	Hutchinson	
5,456,654 A		10/1995	Ball	
5,554,096 A		9/1996	Ball	
5,630,835 A	*	5/1997	Brownlee	607/60
5,716,407 A	*	2/1998	Knapp et al.	128/898
5,749,912 A		5/1998	Zhang	
6,208,235 B1	*	3/2001	Trontelj	340/10.1

OTHER PUBLICATIONS

Brackmann et al., "Evaluation of MRI Compatibility of the
Modified Nucleus Multichannel Auditory Brainstem and

Cochlear Implants," *The American J. of Otology*
17(5):724-9, Sep. 1996.

Teissl et al., "Cochlear Implants: In Vitro Investigation of
Electromagnetic Interference at MR Imaging—Compatibil-
ity and Safety Aspects," *Radiology* 208(3):700-8, Sep.
1998.

Teissl et al., "Magnetic Resonance Imaging and Cochlear
Implants: Compatibility and Safety Aspects," *J. Magn.
Reson. Imaging* 9(1):26-38, Jan., 1999.

* cited by examiner

Primary Examiner—David H. Willse

(74) *Attorney, Agent, or Firm*—Bromberg & Sunstein LLP

(57) **ABSTRACT**

Interference-free coil systems are coil systems having at
least two coils which are identical in terms of inductance.
These coils are arranged such that their magnetic fields are
antiparallel to one another. Consequently, induced voltages
within the coils are substantially eliminated when the coils
are exposed to a homogeneous electromagnetic field. If
exposed to a nonhomogeneous electromagnetic field,
however, a net voltage is induced and enables the extraction
of data and power. Reed switch configurations in the
implantable prostheses protect against induced voltages
caused by the radio frequency field generated by an MR
imager when the reed switches are mounted parallel to the
plane of a receiver. Reed switch configuration may be used
to disconnect, de-tune, or short circuit a receiver. For
example, they may be used to disconnect the receiver
diodes. Some magnet configurations reduce torque caused
by an external magnetic field and prevent demagnetization
when disposed within, outside, or partially within an
implantable prosthesis. Magnets which align with the exter-
nal magnetic field also reduce the torque caused by the field
and prevent demagnetization.

5 Claims, 6 Drawing Sheets

